



866-414-9700

www.TheCPAPShop.com



DeVilbiss Traveler Portable
Compressor Nebulizer Unit
with battery attached
Model 6910P-DR

DeVilbiss Traveler® Portable Compressor Nebulizer System

The **perfect** travel companion

The **DeVilbiss Traveler Portable Compressor Nebulizer System** is designed to help active customers get more out of their aerosol therapy, and **more out of life**. It's never been easier or more convenient to nebulize outside of the home.

The DeVilbiss Traveler includes compressor unit, **disposable** and reusable nebulizers, **five extra filters**, tubing adapter, adult and **pediatric** masks, AC adapter, DC adapter, carrying case, instruction guide and **instructional DVD**, as well as a **full-year** battery warranty. Units can be ordered with **optional battery**.

Performance Features

- Delivers a **consistent and concentrated treatment** of small particles that penetrate deeply into the lungs, delivering medication where it's needed most
- **Fast treatment time** (approximately 10 minutes) with DeVilbiss 900 Series Reusable Nebulizer encourages patient compliance
- Can be used with **any manufacturer's nebulizer** – treatment times will vary
- **45-minute battery life** far outlasts treatment time

Design Features

- Easy, **one-touch operation**
- **Compact and lightweight**
- AC adapter and DC car adapter
- One-year battery and three-year compressor warranty

Model 6910D-DR Includes:

- Compressor unit
- Reusable nebulizer with tubing
- Disposable nebulizer
- 5 replacement air filters and filter door
- Tubing adapter
- Adult and **pediatric masks**
- AC adapter
- DC car adapter
- Easy tote carrying case
- Instructional manual and DVD

Model 6910P-DR Includes:

- All of the above items, plus a rechargeable battery pack

Free Consumer Marketing Program to Generate Sales

- **Comprehensive Traveler marketing kit** includes counter card with reminder tear pad, product flyer (staple to prescription bags) and patient postcard mailers to **generate business**



DeVilbiss Traveler
Portable Compressor
Nebulizer System
Model 6910P-DR

Your nebulizer may be different than shown.

Accessories

6910D-601	Rechargeable Battery
6910D-602	AC Adapter
6910D-603	DC Car Adapter
6910D-605	Air-Inlet Filter (5), Filter Cover (1), & Neb Bracket (1)
6910D-607	Carrying Case
900D-610	Reusable Nebulizer (10/case)
3655D-621	Disposable Nebulizer (50/case)
932168	Instructional DVD
932167	Instructional VHS Video (NTSC format)
AAM12	Adult Aerosol Mask (12/case)
PAM10	Pediatric Aerosol Mask (50/case)
DL1050	Pediatric Character Mask (50/case)



Technical Specifications

MMAD	5 microns	Sound Level	54 dBA
Size w/o Battery	3.34" W x 2.2" H x 3.8" D	Max. Pressure	30 psig or greater
Size w/ Battery	3.34" W x 2.2" H x 5.3" D	Free Air Flow	6 lpm or greater
Weight w/o Battery	.8 lbs	Intermittent Use	30 minutes on/30 minutes off
Weight w/ Battery	1.5 lbs	Master Carton Size	8 units
Electrical	100-240V, 50-60Hz, 12VDC, 1.5A	Battery Life	45 minutes



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com